










	Ameda Finesse™ with Dottie Tote	Ameda® Finesse™ with Minnie Tote	Ameda® Finesse™	Evenflo Deluxe Advanced	Freemie Independence	Lansinoh Smartpump	Lansinoh Signature Pro	Medela Pump in Style Advance® Starter Set	Spectra S2 Plus
									
Unique Features	Finesse mimics the best-in-class Ameda Platinum multi-user hospital pump waveform in a personal pump. This waveform technology helps deliver a smoother, more consistent feel in the vacuum pattern for mom's comfort, which means more milk. Finesse is Ameda's quietest personal electric pump.	Finesse mimics the best-in-class Ameda Platinum multi-user hospital pump waveform in a personal pump. This waveform technology helps deliver a smoother, more consistent feel in the vacuum pattern for mom's comfort, which means more milk. Finesse is Ameda's quietest personal electric pump.	Finesse mimics the best-in-class Ameda Platinum multi-user hospital pump waveform in a personal pump. This waveform technology helps deliver a smoother, more consistent feel in the vacuum pattern for mom's comfort, which means more milk. Finesse is Ameda's quietest personal electric pump.	PerfectlyPure™ design is a closed system that prevents milk back up in tubing, which helps to keep tubing clean and dry, protecting mom's breast milk from contamination and making pump cleaning more efficient. AdvancedControl™ technology creates 32 unique setting combinations of speed and suction for optimal control.	Hospital Grade power that fits in your pocket. The Patented Freemie Independence Pump and new Freemie Closed System Collection Cups (no milk in the tubing!) allows you to double pump, hands free, with your clothes on while you work or play. Save time and gain flexibility by pumping where you need to pump, when you need to pump.	A Smarter Way to Pump! Double Electric Breastpump uses Bluetooth® technology to connect seamlessly to the Lansinoh® Baby App to track pumping sessions, baby's activities and more. Three customizable pumping styles mimic baby's natural feeding pattern. Hygienic closed system.	The Lansinoh® Signature Pro™ Double Electric Breast Pump offers independent control of three different pumping styles and eight levels of suction to maximize milk production and comfort. The closed system design guarantees milk will not back up into the motor and prevents mold and bacteria. It's compact size and quiet motor make it convenient for the mom on-the-go	From the #1 recommended breast pump brand delivers effective pumping in a modular solution. The award-winning Pump in Style has been the first choice of millions of moms for more than two decades. Offering trusted, research based performance with 2-Phase Expression® Technology in a streamlined format customizable to mom's changing lifestyle.	Hospital-strength double breast pump. Hygienic closed system designed to initiate, increase and sustain a milk supply over a long period of time. Gentle suckling action and massage mode to encourage speedy letdown. Adjustable suction strength and cycle speed, night light with two settings, timer
Breast Flanges	2 - 25.0mm, 2 - 30.5mm and 2 - 28.5 inserts	2 - 25mm	2 - 25mm	2 - 30.5mm 2 - 24.5mm and 2 - 28mm Inserts	2 - 25mm (2 - 28mm included Deluxe version only)	2 - 25mm	2 - 25mm	2 - 24mm 2 - 27mm	2 - 24mm 2 - 28mm
Carrying Bag Included	Dottie Tote Bag	Minnie Tote Bag	No	Shoulder Bag	Compact Bag	No	No	No	No
Open or Closed System	Closed	Closed	Closed	Closed	Closed	Closed	Closed	Open	Closed
BPA Free Components	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Bottles Included	4 - 4oz.	2 - 4oz.	4 - 4oz.	4 - 5oz.	None	2 - 5oz	2 - 5oz	2 - 5 oz.	2 - 6 oz.
Other Accessories Included	Dual HygieniKit® Milk Collection System, (2) Extra Valves, Ac Power Adapter, Cool'N Carry Tote with (3) Cooling Elements, Milk Storage Guidelines Card, NoShow Premium Disposable Nursing Pads Sample Pack, Store'N Pour Milk Storage Bags Sample Pack & Instructions for Use	Dual HygieniKit Milk Collection System with (2) 25.0mm CustomFit™ Flanges , (2) Extra Valves, AC Power Adapter, Instructions for Use	Dual HygieniKit Milk Collection System with (2) 25.0mm CustomFit™ Flanges , (2) Extra Valves, AC Power Adapter, Instructions for Use	Insulated Cooler Bag, (3) Ice Packs, Breastfeeding Guide Card, Instructions to Access Simply Breastfeeding™ and Breast Pumps & Briefcases™ Programs , AC Adapter, (2) Extra Valves and Membranes, Instructions for Use	3 Memory Settings, Removable Belt Clip, Wall Adapter	2 Diaphragms and Caps, 2 Breast Flange Bodies, 2 Storage Caps, 2 Feeding Collars, 2 Bottle stands, 2 white valves (2 extra), 3 Connection tubes and tubing connectors, BPA and BPS Free, Pump manual	2 Diaphragms and Caps, 2 Breast Flange Bodies, 2 Storage Caps, 2 Feeding Collars, 2 Bottle stands, 2 white valves (2 extra), 3 Connection tubes and tubing connectors, BPA and BPS Free, Pump manual	(2) Breast Shield Connectors, (2) Valves and membranes, (1) Complete Tubing, (8) Disposable Nursing Pad samples (not pictured), (1) Power adaptor, (1) Battery pack (8 - AA batteries not included) and (1) Instructions for use.	2 valves, 2 sets of tubing and 2 backflow protectors
Battery Compartment	Yes (batteries not included)	Yes (batteries not included)	Yes (batteries not included)	Yes (batteries not included)	Yes - Battery Powered, micro USB Rechargeable (Includes Micro USB Cable)	Yes (batteries not included)	Yes (batteries not included)	Yes (batteries not included)	No
Warranty	2 Year	2 Year	2 Year	1 Year	1 Year	1 Year	1 Year	1 Year	2 Year
Pump Weight	1.1 lbs.	1.1 lbs.	1.1 lbs.	3.15 lbs.	8 oz.	3 lbs.	3 lbs.	4.7 lbs.	4.1 lbs.

The type of breast pump available to you depends upon your insurance coverage. We are contracted with most insurance plans.