








| | Ameda Finesse™ with Dottie Tote | Ameda® Finesse™ with Minnie Tote | Evenflo Deluxe Advanced | Freemie Independence | Lansinoh Smartpump | Medela Pump In Style | Spectra S2 Plus |
|-----------------------------------|---|---|--|---|---|---|--|
| |  |  |  |  |  |  |  |
| Unique Features | Finesse mimics the best-in-class Ameda Platinum multi-user hospital pump waveform in a personal pump. This waveform technology helps deliver a smoother, more consistent feel in the vacuum pattern for mom's comfort, which means more milk. Finesse is Ameda's quietest personal electric pump. | Finesse mimics the best-in-class Ameda Platinum multi-user hospital pump waveform in a personal pump. This waveform technology helps deliver a smoother, more consistent feel in the vacuum pattern for mom's comfort, which means more milk. Finesse is Ameda's quietest personal electric pump. | PerfectlyPure™ design is a closed system that prevents milk back up in tubing, which helps to keep tubing clean and dry, protecting mom's breast milk from contamination and making pump cleaning more efficient. AdvancedControl™ technology creates 32 unique setting combinations of speed and suction for optimal control. | Hospital Grade power that fits in your pocket. The Patented Freemie Independence Pump and new Freemie Closed System Collection Cups (no milk in the tubing!) allows you to double pump, hands free, with your clothes on while you work or play. Save time and gain flexibility by pumping where you need to pump, when you need to pump. | A Smarter Way to Pump! Double Electric Breastpump uses Bluetooth® technology to connect seamlessly to the Lansinoh® Baby App to track pumping sessions, baby's activities and more. Three customizable pumping styles mimic baby's natural feeding pattern. Hygienic closed system. | Patented 2-Phase Expression technology pump. Adjustable speed/vacuum control helps mom choose her most comfortable pump setting. One-touch let-down button lets mom toggle between Stimulation and Expression Phases for more efficient pumping sessions. | Hospital-strength double breast pump. Hygienic closed system designed to initiate, increase and sustain a milk supply over a long period of time. Gentle suckling action and massage mode to encourage speedy letdown. Adjustable suction strength and cycle speed, night light with two settings, timer |
| Breast Flanges | 2 - 25.0mm, 2 - 30.5mm and 2 - 28.5 inserts | 2 - 25mm | 2 - 30.5mm 2 - 24.5mm and 2 - 28mm Inserts | 2 - 25mm (2 - 28mm included Deluxe version only) | 2 - 25mm | 2 - 24 mm | 2 - 24mm and 2 - 28mm |
| Carrying Bag Included | Dottie Tote Bag | Minnie Tote Bag | Shoulder Bag | Compact Bag | No | No | No |
| Open or Closed System | Closed | Closed | Closed | Closed | Closed | Open | Closed |
| BPA Free Components | Yes | Yes | Yes | Yes | Yes | Yes | Yes |
| Bottles Included | 4 - 4oz. | 2 - 4oz. | 4 - 5oz. | None | 2 - 5oz | 2 - 5 oz. | 2 - 6 oz. |
| Other Accessories Included | Dual HygieniKit® Milk Collection System, (2) Extra Valves, Ac Power Adapter, Cool'N Carry Tote with (3) Cooling Elements, Milk Storage Guidelines Card, NoShow Premium Disposable Nursing Pads Sample Pack, Store'N Pour Milk Storage Bags Sample Pack & Instructions for Use | Dual HygieniKit Milk Collection System, (2) Extra Valves, AC Power Adapter & Instructions for Use | Insulated Cooler Bag, (3) Ice Packs, Breastfeeding Guide Card, Instructions to Access Simply Breastfeeding™ and Breast Pumps & Briefcases™ Programs, AC Adapter, (2) Extra Valves and Membranes, Instructions for Use | 3 Memory Settings, Removable Belt Clip, Wall Adapter | 2 Diaphragms and Caps, 2 Breast Flange Bodies, 2 Storage Caps, 2 Feeding Collars, 2 Bottle stands, 2 white valves (2 extra), 3 Connection tubes and tubing connectors, BPA and BPS Free, Pump manual | Compact motor unit in a soft bag. | 2 valves, 2 sets of tubing and 2 backflow protectors |
| Battery Compartment | Yes (batteries not included) | Yes (batteries not included) | Yes (batteries not included) | Yes - Battery Powered, micro USB Rechargeable (Includes Micro USB Cable) | Yes (batteries not included) | No | No |
| Warranty | 2 Year | 2 Year | 1 Year | 1 Year | 1 Year | 1 Year | 2 Year |
| Pump Weight | 1.1 lbs. | 1.1 lbs. | 3.15 lbs. | 8 oz. | 3 lbs. | 4.2 lbs. | 4.1 lbs. |

The type of breast pump available to you depends upon your insurance coverage. We are contracted with most insurance plans.