







| | Ameda Purely Yours Ultra | Evenflo Deluxe Advanced | Freemie Freedom | Lansinoh Smartpump | Medela Pump In Style | Spectra S2 Plus |
|-----------------------------------|---|--|---|---|---|--|
| |  |  |  |  |  |  |
| Unique Features | Proven Airlock Protection™ barrier which helps protect milk while keeping pump tubing dry and clean. No need to clean the tubing. CustomControl™ empowers moms to more closely mimic their baby's nursing rhythm. Moms have complete control of their cycle, both suction and speed, to select the best levels for their milk flow and comfort. | PerfectlyPure™ design is a closed system that prevents milk back up in tubing, which helps to keep tubing clean and dry, protecting mom's breast milk from contamination and making pump cleaning more efficient. AdvancedControl™ technology creates 32 unique setting combinations of speed and suction for optimal control. | Allows you to double pump hands free and with your shirt on. The Freemie cups are supported by regular or nursing bras. Pump connection kit allows you to interrupt pumping and disconnect from pump so you can move around without interrupting milk collection, and then reconnect and resume pumping | A Smarter Way to Pump! Double Electric Breastpump uses Bluetooth® technology to connect seamlessly to the Lansinoh® Baby App to track pumping sessions, baby's activities and more. Three customizable pumping styles mimic baby's natural feeding pattern. Hygienic closed system. | Patented 2-Phase Expression technology pump. Adjustable speed/vacuum control helps mom choose her most comfortable pump setting. One-touch let-down button lets mom toggle between Stimulation and Expression Phases for more efficient pumping sessions. | Hospital-strength double breast pump. Hygienic closed system designed to initiate, increase and sustain a milk supply over a long period of time. Gentle suckling action and massage mode to encourage speedy letdown. Adjustable suction strength and cycle speed, night light with two settings, timer |
| Breast Flanges | 2 - 25mm and 2 - 30.5mm and 2 - 28.5mm inserts | 2 - 30.5mm 2 - 24.5mm and 2 - 28mm Inserts | 2 - 25mm 2 - 28mm (Deluxe version only) | 2 - 25mm | 2 - 24 mm | 2 - 24mm and 2 - 28mm |
| Carrying Bag Included | Dottie Tote | Shoulder Bag | No | No | No | No |
| Open or Closed System | Closed | Closed | Closed | Closed | Open | Closed |
| BPA Free Components | Yes | Yes | Yes | Yes | Yes | Yes |
| Bottles Included | 4 - 4 oz. | 4 - 5 oz. | None | 2 - 5 oz. | 2 - 5 oz. | 2 - 6 oz. |
| Other Accessories Included | Cool'N Carry Tote with (3) cooling elements, (4) extra 4 oz. bottles with lock-tight sealing lids, Milk Storage Guidelines card, NoShow Premium™ Disposable Nursing Pads 2 pack Sample, Store'N Pour™ Milk Storage Bags 2 pack Sample, (2) extra valves | Insulated Cooler Bag, (3) Ice Packs, Breastfeeding Guide Card, Instructions to Access Simply Breastfeeding™ and Breast Pumps & Briefcases™ Programs, AC Adapter, (2) Extra Valves and Membranes, Instructions for Use | Integrated powercord | 2 Diaphragms and Caps, 2 Breast Flange Bodies, 2 Storage Caps, 2 Feeding Collars, 2 Bottle stands, 2 white valves (2 extra), 3 Connection tubes and tubing connectors, BPA and BPS Free, Pump manual | Compact motor unit in a soft bag. | 2 valves, 2 sets of tubing and 2 backflow protectors. |
| Battery Compartment | Yes (6-AA batteries not included) | Yes (6-AA batteries not included) | No | Yes (6-AA batteries not included) | No | No |
| Warranty | 1 Year | 1 Year | 1 Year | 1 Year | 1 Year | 2 Year |
| Pump Weight | 4.0 lbs | 3.15 lbs | 3.2 lbs | 3 lbs | 4.2 lbs | 4.1 lbs |

The type of breast pump available to you depends upon your insurance coverage. We are contracted with most insurance plans.