








	Ameda Finesse™ with Dottie Tote	Ameda® Finesse™ with Minnie Tote	Evenflo Deluxe Advanced	Freemie Freedom	Lansinoh Smartpump	Medela Pump In Style	Spectra S2 Plus
							
Unique Features	Finesse mimics the best-in-class Ameda Platinum multi-user hospital pump waveform in a personal pump. This waveform technology helps deliver a smoother, more consistent feel in the vacuum pattern for mom's comfort, which means more milk. Finesse is Ameda's quietest personal electric pump.	Finesse mimics the best-in-class Ameda Platinum multi-user hospital pump waveform in a personal pump. This waveform technology helps deliver a smoother, more consistent feel in the vacuum pattern for mom's comfort, which means more milk. Finesse is Ameda's quietest personal electric pump.	PerfectlyPure™ design is a closed system that prevents milk back up in tubing, which helps to keep tubing clean and dry, protecting mom's breast milk from contamination and making pump cleaning more efficient. AdvancedControl™ technology creates 32 unique setting combinations of speed and suction for optimal control.	Allows you to double pump hands free and with your shirt on. The Freemie cups are supported by regular or nursing bras. Pump connection kit allows you to interrupt pumping and disconnect from pump so you can move around without interrupting milk collection, and then reconnect and resume pumping	A Smarter Way to Pump! Double Electric Breastpump uses Bluetooth® technology to connect seamlessly to the Lansinoh® Baby App to track pumping sessions, baby's activities and more. Three customizable pumping styles mimic baby's natural feeding pattern. Hygienic closed system.	Patented 2-Phase Expression technology pump. Adjustable speed/vacuum control helps mom choose her most comfortable pump setting. One-touch let-down button lets mom toggle between Stimulation and Expression Phases for more efficient pumping sessions.	Hospital-strength double breast pump. Hygienic closed system designed to initiate, increase and sustain a milk supply over a long period of time. Gentle suckling action and massage mode to encourage speedy letdown. Adjustable suction strength and cycle speed, night light with two settings, timer
Breast Flanges	2 - 25.0mm, 2 - 30.5mm and 2 - 28.5 inserts	2 - 25mm	2 - 30.5mm 2 - 24.5mm and 2 - 28mm Inserts	2 - 25mm 2 - 28mm (Deluxe version only)	2 - 25mm	2 - 24 mm	2 - 24mm and 2 - 28mm
Carrying Bag Included	Dottie Tote Bag	Minnie Tote Bag	Shoulder Bag	No	No	No	No
Open or Closed System	Closed	Closed	Closed	Closed	Closed	Open	Closed
BPA Free Components	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Bottles Included	4 - 4oz.	2 - 4oz.	4 - 5 oz.	None	2 - 5 oz.	2 - 5 oz.	2 - 6 oz.
Other Accessories Included	Dual HygieniKit® Milk Collection System, (2) Extra Valves, AC Power Adapter, Cool'N Carry Tote with (3) Cooling Elements, Milk Storage Guidelines Card, NoShow Premium Disposable Nursing Pads Sample Pack, Store'N Pour Milk Storage Bags Sample Pack & Instructions for Use	Dual HygieniKit Milk Collection System, (2) Extra Valves, AC Power Adapter & Instructions for Use	Insulated Cooler Bag, (3) Ice Packs, Breastfeeding Guide Card, Instructions to Access Simply Breastfeeding™ and Breast Pumps & Briefcases™ Programs , AC Adapter, (2) Extra Valves and Membranes, Instructions for Use	Integrated powercord	2 Diaphragms and Caps, 2 Breast Flange Bodies, 2 Storage Caps, 2 Feeding Collars, 2 Bottle stands, 2 white valves (2 extra), 3 Connection tubes and tubing connectors, BPA and BPS Free, Pump manual	Compact motor unit in a soft bag.	2 valves, 2 sets of tubing and 2 backflow protectors.
Battery Compartment	Yes (6-AA batteries not included)	Yes (6-AA batteries not included)	Yes (6-AA batteries not included)	No	Yes (6-AA batteries not included)	No	No
Warranty	2 Year	2 Year	1 Year	1 Year	1 Year	1 Year	2 Year
Pump Weight	1.1 lbs.	1.1 lbs.	3.15 lbs	3.2 lbs	3 lbs	4.2 lbs	4.1 lbs

The type of breast pump available to you depends upon your insurance coverage. We are contracted with most insurance plans.